C00P1484_V08

BLUE MARKER AGUETTANT®

10 mg/mL

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Fabricant

Laboratoire AGUETTANT
1, rue Alexander Fleming - 69007 LYON - France

Dispositif

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BLUE MARKER AGUETTANT® 10 mg/mL
Methylene blue

Indication

Visualisation aid for surgical procedures:

- Demarcation of surgical tissues and operative specimens
- Seal test for urinary and colorectal sutures
- Evaluation of tubal patency

Target treated population

BLUE MARKER AGUETTANT® is intended to be used in paediatric and adult populations (excluding pregnant and breastfeeding women), who do not meet contraindications, and regardless of gender.

Intended users

BLUE MARKER AGUETTANT® is intended to be used by healthcare professionals.

Composition

Methylene blue, 3H₂O (Methylthionium chloride, 3H₂O) 1 g
Water for injections *q.s.* 100 mL
Sodium hydroxide and/or hydrochloric acid *q.s.* pH 3.0-4.5

Clinical benefit / Performance / Mode of actions

Steam sterilization

Clinical benefit

- BLUE MARKER AGUETTANT® allows to limit excision to tissue marked with dye,
- BLUE MARKER AGUETTANT® allows routine per-operative visualisation of any leak from urinary and colorectal sutures,
- BLUE MARKER AGUETTANT® allows to confirm tubal patency via confirmation of the flow of diluted methylene blue through the Fallopian tubes in the context of woman infertility.

Device Performance

BLUE MARKER AGUETTANT® is a visualisation aid for surgical procedures:

- Demarcation of surgical tissues and operative specimens:
 - Identification of the renal cavity during percutaneous renal surgery (nephrolithotomy)
 - Intraoperative identification of a urinary or colonic fistula before excision
 - Intraoperative demarcation of the pilonidal sinus before removal
 - Marking the future opening of a stoma before excision

- Seal test for urinary and colorectal sutures
- Evaluation of tubal patency

Requirement before Use and Instructions for use**Method of administration and dosage**

For all indications, BLUE MARKER AGUETTANT® should be diluted 1 part to 100 parts before use with a sterile water solution or in 0.9% sodium chloride sterile solution (to produce a solution of 0.01% methylene blue). The solution should be used immediately after dilution.

The method of administration and the volume of diluted solution to be administered vary depending on the nature and size of the area to be visualised.

The following volumes are for information only.

Demarcation of surgical tissues and operative specimens:

- Identification of the renal cavity during percutaneous renal surgery (nephrolithotomy): Intravesically administer the diluted solution (100 ml to 250 ml commonly used).
- Intraoperative identification of a urinary or colonic fistula before excision: Intravesically administer the diluted solution (10 ml to 100 ml commonly used).
- Intraoperative demarcation of the pilonidal sinus before removal: Use a syringe to administer the diluted solution directly into the opening of the pilonidal sinus (volume is based on the size of the cyst, generally 1 ml to 5 ml).
- Marking the future opening of a stoma before excision: The volume of diluted solution to be administered via local subcutaneous injection is at the discretion of the hospital doctor (2 ml to 5 ml commonly used).

Seal test for sutures during surgery:

- Urinary and colorectal sutures: Intravesically or rectally administer the diluted solution (25 ml to 250 ml depending on the location of the suture). Lack of blue colour at the suture zone validates the suture.

Evaluation of tubal patency:

- Administer the diluted solution (10 ml commonly used) via the intrauterine route. The flow of diluted methylene blue through the Fallopian tubes confirms tubal patency.

The placement of white sterile compresses in contact with the area to be tested, depending on the indications, can assist in the visualisation of a leak of diluted methylene blue

Warnings, precautions for use and contraindications**Contraindication**

Do not use BLUE MARKER AGUETTANT®:

- in case of treatment with drug with serotonergic activity including selective serotonin reuptake inhibitors (SSRIs), serotonin-norepinephrine reuptake inhibitors (SNRIs: venlafaxine, desvenlafaxine, duloxetine, fluoxetine), tricyclic antidepressants (clomipramine, amitriptyline), and other drugs with serotonergic activity such as bupropion, buspirone, trazodone, fentanyl, meperidine, tramadol, methadone. Use with caution in case of concomitant treatment with others Monoamine Oxidase Inhibitors.
- in case of glucose-6-phosphate dehydrogenase deficiency,
- in case of known hypersensitivity to methylene blue or other thiazine dyes,
- in case of pregnancy and breastfeeding.
- Use with caution BLUE MARKER AGUETTANT®:
- in case of moderate or severe kidney disease.

Precaution for use

When possible, rinse the stained area with the same solvent used for the dilution. BLUE MARKER AGUETTANT® must be used by healthcare professionals. It is recommended to wear gloves during use.

Do not use if the bottle or the cap have been damaged.

Single use. Risk of microbial contamination for the patient in the event of re-use.

Do not use after the date which is stated on the carton.

Warning

Do not inject BLUE MARKER AGUETTANT® intravenously, intrathecally, intra-amniotically, intraocularly and orally.

Avoid contact with the surface of the eye.

In case of contact with the ocular surface, rinse thoroughly with water.

Do not use if crystallization reaction has occurred in the diluted methylene blue solution and if the solution is colourless.

Undesirable effects

Any serious incident related to the device should be notified to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established. Undesirable effects vary depending on the use of the medical device. The following effects have been reported:

- hypersensitivity reaction, anaphylactic reaction, pulmonary oedema,
- intraoperative - or postoperative methemoglobinemia,
- vasoconstriction, ulceration, necrosis,
- temporary discolouration of the urine, faeces, skin, hair and nails.
- False positive or false negative, in the context of seal test for urinary sutures. The use of methylene blue in patients with glucose-6-phosphate dehydrogenase deficiency may accelerate haemolysis.
- The use of methylene blue in combination with selective serotonin reuptake inhibitors (SSRIs) may cause severe central nervous system toxicity. Subcutaneous injection may cause vasoconstriction or skin necrosis.

Storage, handling and discard

There is no special storage or handling conditions. Do not refrigerate. Do not Freeze. Dispose of the vial according to current practices.

Information if the patient is not the user

The Healthcare Professional should inform the patient regarding the undesirable effects and the contraindications.

Symbols

 Sterilized by steam sterilization
 Batch code
 Use-by date
 Single-use
 Consult instructions for use
 Manufacturer

Date of revision of the text:

Version 02 – 11/2023

CE 0459
(2011)

AGUETTANT