

Innovation for all

AGUETTANT® SELF-FLUSHING INFUSION BAG *PLUS*



| A NEW-GENERATION OF INFUSION BAG FOR HOSPITAL AND HOME CARE |

AGUETTANT® SELF-FLUSHING INFUSION BAG PLUS

A revolutionary concept

With its dual flushing, Aguettant's infusion bag **PLUS** enables to purge the line prior to the drug infusion and to perform an automatic flushing of the connecting line once the drug is infused without any intervention from the health care professional.


Instructions of use:



1 - Self-Flushing Infusion bag before use



2 - Purge the line



3 - Break the seal



4 - Open the line to infuse the drug



5 - Bag flushes automatically once the infusion is completed

worldwide
patent pending



Solutions for preserving life

Features of Self-flushing infusion bag PLUS:

→ Full dose delivery : an improved compliance

A study* demonstrated that the use of conventional infusion bags induces an under-dosage and therapy inefficacy risks, as up to 20% of the drug is lost in the connecting line after the drug infusion. With Aguettant® Self-Flushing infusion bag, the full dose is systematically delivered without any intervention from the medical staff. (*ANFAR 27 52008) 514-519)

→ Maximum safety : a reduced risk of accidents and infections

PLUS: The purge of the line prior to the delivery of the drug, and the automatic post flushing eliminate the risk of exposure to toxic drug during the line connection and disconnection. The system is totally closed.

→ Fast and convenient use

No change of the current infusion protocol.
No need for the nurse to intervene at the end of the infusion.

Benefits:

→ For the patient:

- Guarantee of compliance to the medical prescription
- Reduced risks of contamination due to minimized manual intervention

→ For health care professionals:

- Time saved
- Maximized safety against risk of exposure to toxic drugs
- Simple and ergonomic design – easy to use

→ For the environment:

- Reduced waste by eliminating the need for syringes, vials of solution and sterile compresses traditionally required
- No drug waste

→ Pharmaceutical industry:

- Patented internationally
- Mechanically reliable and safe
- Adaptable to different infusion bags size
- Dual chamber infusion bag also available
- Exclusive licence available

www.aguettant.com

Laboratoire **AGUETTANT**
Headquarters & Lyon manufacturing plant
1 Rue Alexander Fleming
69007 Lyon